

TOMEY TABLES

ELECTRIC LIFT TABLES

DELIGHT IN SIGHT

Shapely design. For everybody.



TTUD-1000



TT-4060

- Individual solution for all products
- Easy to install
- Silent & smoothly working lifting column
- Stable & comfortable
- Shelves available for cables & accessories
- Equipped with 4 casters

 **TOMEY**
TECHNOLOGY AND VISION

TOMEY ELECTRIC LIFT TABLES



QUALITY IN DETAIL

Find your perfect solution for all products – all individual demands can be fulfilled with one of our table series. If ultrasound, refraction instruments or a high class OCT, our **TOMEY tables** are built to offer a perfect match to your needs. All our models are slimline designed, equipped with a very robust base, known as one of the most stable, silent and smoothly working ophthalmic tables on the market. The surface of the table tops is easy to clean as this is a must in the medical field. Provided with four casters made to be used under pressure, each table can be easily moved with the instruments on top. Clever solution for a smooth operation!

BASIC LINE



Three vertically adjustable electric lift tables that allow the use of instruments up to a weight of 80 kg. Optional shelf TTPR-1000 for OA-2000 printer (only for TT-1060) is available.

TT2C-1000



TT2C-1000 is a vertically adjustable electric table with two columns. This makes it ideal for operation with two systems (e.g. FT+RC or TMS+OA). Additionally it is very useful for wheel chair access. The maximum load is 150 kg. Optional shelves for cables and accessories are available.

TTUD-1000



This extremely loadable model for our ultrasound line features one drawer (for storage of the probes), holder for tissues and cleaning fluids. The maximum load is 80 kg. Optimised space for our recommended Mitsubishi P95DE video printer (USB).

SPECIFICATIONS

PRODUCT	TABLE TOP DIMENSIONS	TABLE TOP HEIGHT
TT-4060	720 x 460 mm	725 - 925 mm
TT-1060	1060 x 600 mm	730 - 930 mm
TTVS-1000	1040 x 553 mm V-shape	725 - 925 mm
TT2C-1000	1240 x 550 mm	609 - 909 mm
TTUD-1000	580 x 545 mm	705 - 1005 mm
OPTIONAL SHELVES		
TTPR-1000 Shelf for OA-2000 printer (only for TT-1060)		
TTCS-1000 Shelf of CASIA adapted shelves (only for TT2C-1000)		
TTWS-1000 Shelf for cable, adapter, connector (only for TT2C-1000)		

2017/08 - subject to change without notice



TOMEY EUROPE
TOMEY GmbH
 Wiesbadener Straße 21
 90427 Nürnberg, Germany
 Phone +49 911 938 546 2 0
 Fax +49 911 938 546 2 20
 Email info@tomey.de

TOMEY ASIA-PACIFIC
TOMEY CORPORATION JAPAN
 2-11-33 Noritakeshinmachi
 Nishi-ku, Nagoya 451-0051, Japan
 Phone +81 52 581 5327
 Fax +81 52 561 4735
 Email intl@tomey.co.jp

 **TOMEY**
 TECHNOLOGY AND VISION

www.tomey.de