

SPECTRA ULTRA
TECHNICAL CHARACTERISTICS

- › Maximum reflection of 0.2 % over the full spectrum of visible light.
- › Maximum transmission 99.8 %.
- › 50 % lower residual color intensity of anti-reflective treatment.
- › Reduced reflection of UV light from the back-side into the eyes.

ANTI-REFLECTIVE TREATMENT

Spectra Ultra	MAX Ultra	Satin	Oasis
 Maximum glare reduction	 High glare reduction	 High glare reduction	 Standard glare reduction
 Maximum resistance to scratches	 High resistance to scratches	 High resistance to scratches	 Standard Resistance to scratches
 UV reduction from the back side	 UV reduction from the back side		
 Premium superhydrophobic layer	 Superhydrophobic layer	 Superhydrophobic layer	 Hydrophobic layer
 Antistatic layer	 Antistatic layer	 Antistatic layer	
 Premium oleophobic layer	 Standard oleophobic layer	 Standard oleophobic layer	
 Residual color: blue/purple	 Residual color: blue/green	 Residual color: green	 Residual color: green

Index 1.5 does not block UV, in case of 100% protection it is necessary to select index 1.6 and higher, or 1.5 UV+ 420 BlueCUT or order a UV filter.



spectra
ultra

PURE BEAUTY.
RESISTANT
AT FIRST
SIGHT.

Thanks to an innovative anti-reflective coating, Spectra Ultra lenses are twice as scratch resistant. They're almost completely clear.


NEW PREMIUM ANTI-REFLECTIVE TREATMENT

Thanks to the innovative anti-reflective coating, **Spectra Ultra** lenses are **twice as scratch resistant** compared to conventional AR treatments.

Spectra Ultra features a **low reflection of 0.2 %** over the full spectrum of visible light.

This feature ensures high contrast vision, glare reduction (in day and night vision) and high lens aesthetics with a very light residual blue-violet color.

SPECTRA VS. SPECTRA ULTRA



*The graph shows the adjustment of the front side of the lens.

BENEFITS

- › The clearest premium quality anti-reflective technology in the Omega Optix portfolio.
- › Low residual color intensity ensuring high aesthetics.
- › Maximum hardness to ensure the scratch resistance.
- › Long-lasting clarity, resistant against dust, humidity, smudges and fingerprints.
- › The highest level of hydrophobic and oleophobic top coating for easy cleaning.

FOR WHOM?

The perfect choice for people who appreciate aesthetics and need a reliable solution with exceptional optical features.

SPECTRA ULTRA IS IDEAL FOR CHALLENGES OF EVERYDAY LIFE.

